



XIANNING FULL GUARD MEDICAL PRODUCTS CO.,LTD.

Add : East Yongan Avenue, Wenquan Economic Development Zone, Xianning City, Hubei Province, 437100, China.








Contact person : Mrs Shirley Yang







E-mail: shirley.yang@fullguard-medical.com www.fullguardmedical.com

Tel: +86-715-8200113 Mob:+86-15377205746 WhatsApp:+8615377205746


Full Guard Disposable Surgical Kits

Item Name	Picture	NO.	Products Name	Material	Size(cm)	Unit(pc)	Packing	FOB Price (USD)	MOQ (set)
品名	图片	编号	内置品名	材质	尺寸	数量	包装方式	FOB Shanghai,China	起订量(包)
通用包 Universal Kits FJB0035		1	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190	1	1set/blister bag 10bags/ctn Eo sterile 47*32*58cm,	7.109USD/PC	2000
		2	Head Drape 加强型头单	35g SMS, 240cm side Adhesive 5*80cm, Reinforced 80g 30*90cm	150*240	1			
		3	Foot Drape 加强型脚单	35g SMS, 180cm side Adhesive 5*80cm, Reinforced 80g 30*90cm	175*180	1			
		4	Side Drape 加强型边单	35g SMS, 90cm side Adhesive 5*88cm, Reinforced 80g 30*90cm	75*90	2			
		5	Mayo Stand Cover 美亚袋	Width 80cm, 5.5s PE film+ reinforced 30g Hydrophilic PP 60*100cm	80*145cm	1			
		6	Op-tape 背胶料	-	10*50	2			
		7	Hand Towels 擦手纸	Two folded	30*40	4			
		8	Suture bag 帆袋	80g Kraft paper	16*28cm	1			
		9	Standard Surgical Gown 标准型手术衣	40g SMS, Ultrasonic welding, knitted cuff, verco neck, book folding	XL/130*150cm	2			
		10	Wrapper 外包布	35g SMS	100*100cm	1			
通用包 Universal Kits FJB0035-01		1	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190cm	1	1set/blister bag 10bags/ctn Eo sterile 47*32*33cm,	4.496USD/PC	2000
		2	Head Drape 加强型头单	35g SMS, 240cm side Adhesive 5*80cm, Reinforced 80g 30*90cm	150*240cm	1			
		3	Foot Drape 加强型脚单	35g SMS, 180cm side Adhesive 5*80cm, Reinforced 80g 30*90cm	190*190cm	1			
		4	Side Drape 加强型边单	35g SMS, 90cm side Adhesive 5*88cm, Reinforced 80g 30*90cm	75*90cm	2			
		6	Op-tape 背胶料	-	10*50	1			
		7	Hand Towels 擦手纸	Two folded, 60g Spunlace	30*40	2			

顺产包 Delivery Pack-FJB0034-01 	1	Bottom Drape 臀底单	6s PE film, with Pouch	100*112cm	1	1set/blister bag 10bags/ctn Eo sterile 50*33*53cm	6.053USD/PC	2000
	2	Baby cloth 婴儿毯	60g Spunlace Material	75*90cm	1			
	3	Leggings 腿套	Blue 40g SMS	75*120cm	2			
	4	Standard Surgical Gown 标准型手术衣	40g SMS, Ultrasonic welding, knitted cuff, verco neck, book folding	XL/130*150cm	2			
	5	Instrument Table Cover 桌布	Width 150cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*200cm	150*200cm	1			
	6	Medium Drape 中单	Blue 40g SMS, Adhesive 3*80cm	150*180cm	1			
	7	Hand Towels 擦手纸	60g Spunlace Material, White	30*40cm	2			
一次性腹部手术包 Laparotomy Pack 	1	Laparoscopy Drape 主单	Blue 35g SMS, Fenestration 10*30cm with Adhesive 2.5cm, Reinforced pad: 55g Water repellent Film	190*1240cm	1	1set/blister bag 10bags/ctn Eo sterile 56*45*46cm	7.424USD/PC	2000
	2	Drape Sheet 铺单	Blue 35g SMS, No adhesive	140*160cm	2			
	3	Utility Drape 简易单	Blue 35g SMS, adhesive 1.2*58cm	40*60cm	4			
	4	Mayo Stand Cover 美亚袋	Blue 35g SMS, reinforced 58g Film 50*55cm	58*138cm	1			
	5	Hand Towels 擦手纸	53gsm white paper	30*40cm	3			
	6	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190cm	1			
	7	Reinforced Surgical Gown 加强手术衣	35g SMS, Ultrasonic welding, reinforced surgical gown, knitted cuff, verco neck, book folding, reinforced pad: 40g PP+PE film	XL/130*150cm	2			
牙科包 Dental Pack-FJB0028 	1	Standard Surgical Gown 标准型手术衣	40g SMS, Ultrasonic welding, knitted cuff, verco neck, book folding	L/120*150cm	2	1set/blister bag 10bags/ctn Eo sterile 47*25*46cm	5.685USD/PC	2000
	2	Op-tape 背胶料	Folded, fixed by the transparent elastic	5*50cm	2			
	3	Hand Towels 擦手纸	53gsm white paper	30*40cm	2			
	4	Gauze 无纺布片	5*5 * 4 layer	5*5cm	10			
	5	Video Camera Cover PE 保护套	3s transparent PE, with cardboard	8*120cm	2			
	6	Aperture Drape 洞巾	Blue, 60g Lamination, Hole 20*20cm	75*90cm	1			
	7	Drape Sheet 铺单	Blue, 60g Lamination, No Adhesive	50*60cm	2			
	8	"U" Sharpe Dental Drape	Blue, 60g Lamination, U sharpe 10*10cm, adhesive 5*5cm around	100*150cm	1			
	9	Suture bag 帆袋	80g Kraft paper	16*28cm	1			
	10	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190cm	1			
	11	Wrapper 外包布	30g SMS	100*100cm	1			
牙科包 Dental Pack-FJB0029 	1	Standard Surgical Gown 标准型手术衣	40g SMS, Ultrasonic welding, knitted cuff, verco neck, book folding	L/120*150cm	2	1set/blister bag 10bags/ctn Eo sterile 47*25*45cm	5.779USD/PC	2000
	2	Op-tape 背胶料	Folded, fixed by the transparent elastic	5*50cm	2			
	3	Hand Towels 擦手纸	53gsm white paper	30*40cm	4			
	4	Gauze 无纺布片	5*5 * 4 layer	5*5cm	10			
	5	Video Camera Cover PE 保护套	3s transparent PE, with cardboard	8*120cm	2			
	6	Drape Sheet 铺单	Blue, 60g Lamination, No Adhesive	75*90cm	2			
	7	"U" Sharpe Dental Drape	Blue, 60g Lamination, U sharpe 10*10cm, adhesive 5*5cm around	100*150cm	1			
	8	Suture bag 帆袋	80g Kraft paper	16*28cm	1			
	9	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190cm	1			
	10	Wrapper 外包布	30g SMS	100*100cm	1			
眼科包 Ophthalmic pack-FJB0027 	1	Instrument Table Cover 桌布	Width 100cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*150cm	100*150cm	1	1set/blister bag 20bags/ctn Eo sterile 47*25*38cm	3.588USD/PC	2000
	2	Ophthalmic drape 眼科单	40g SMS, Pouch: 25*30cm, Incise: 15*18cm, Hole: 9*12cm	150*200cm	1			
	3	Hand Towels 擦手纸	60g Spunlace Material, White	30*40cm	2			
	4	Bundle Cover 外包布	30g SMS	100*100cm	1			
眼科包 Ophthalmic pack 	1	Instrument Table Cover 桌布	Width 100cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*150cm	100*150cm	1	1set/blister bag 20bags/ctn Eo sterile 47*25*40cm	3.782USD/PC	2000
	2	Single Pouch Ophthalmic drape 眼科单	40g SMS, Middle fenestration 12cm, 3M surgical film 15*15cm+1pc 8S transparent PE Pouch	100*130cm	1			
	3	Standard Surgical Gown 标准型手术衣	40g SMS, Ultrasonic welding, knitted cuff, verco neck, book folding	XL/130*150cm	2			
	4	Hand Towels 擦手纸	60g Spunlace Material, White	30*40cm	2			
	5	Bundle Cover 外包布	30g SMS	100*100cm	1			
血管造影包 Angio Pack-FJB0019 	1	Angio Drape 造影包	40g SMS, 2 Fenestrations	210*300cm	1	1set/blister bag 8bags/ctn Eo sterile 51*45*42cm,	11.55USD/PC	2000
	2	Hand towels 擦手纸	60g Spunlace 60g	30*40cm	2			
	3	Reinforced Surgical Gown 加强型手术衣	40g SMS, Ultrasonic welding, knitted cuff, verco neck, book folding	XL/130*150cm	2			
	4	Circular Banded Bag 圆形仪器挡板套	Transparent PE	diameter100cm	1			
	5	Rectangular Banded Bag 方形仪器挡板套	Transparent PE	75*90cm	1			
	6	Solution Cup 量杯	PP	500cc	1			
	7	X-Ray Gauze	7.5*7.5cm-8ply with X-Ray	7.5*7.5cm	10			
	8	Betadine Spone 海绵棒		ST	1			
	9	Gallipots 量杯	PP	1000cc	1			
	10	Guid Wire 导丝管	0.035in	150cm	1			
	11	Instrument Table Cover 桌布	Width 150cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	150*190cm	1			
	12	SMS Wrap 包布	30g SMS	100*100cm	1			

Cardiovascular Pack 心血管手术包		1	Cardiovascular U Split Drape/主单	40g SMS	215/275x335cm	1	1set/Blister bag, 6sets/Ctn, EO Sterile,Ctn size: 60x40x45cm	15.28USD/PC	2000
		2	Bilateral Split Drape/双侧分体单	55g PE film	200*200cm	1			
		3	Side Drape 加强型边单	40g SMS, 90cm side Adhesive 5*88cm, Reinforced 80g 30*90cm	75*90	1			
		4	Mayo Stand Cover美亚袋	Width 80cm, 5.5s PE film+ reinforced 30g Hydrophilic PP 60*100cm	80*145cm	1			
		5	Hand Towels 擦手纸	Two folded.60g Spunlace material, 60g 水刺布	30*40	4			
		6	Reinforced Surgical Gown 加强手术衣	40g SMS, Ultrasonic welding, reinforced surgical gown, knitted cuff, verco neck, book	XL/130*150cm	2			
		7	Op-tape背胶料	-	10*50	2			
		8	Suture bag 幘袋	80g Kraft paper	16*28cm	2			
		9	Instrument Table Cover 桌布	Width 150cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*230cm	150*230	1			
		10	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190	1			
		11	Wrapper 外包布	30g SMS	100*100cm	1			
ENT Pack/耳鼻喉科手术包		1	Split Drape/分切主单	40g SMS	200*300cm	1	1set/Blister bag, 10sets/Ctn, EO Sterile,Ctn size: 60x40x45cm	7.670USD/PC	2000
		2	Head Bar Drape/头巾	40g SMS	100*100cm	1			
		3	Mayo Stand Cover美亚袋	Width 80cm, 5.5s PE film+ reinforced 30g Hydrophilic PP 60*20cm	80*145cm	1			
		4	Hand Towels 擦手纸	60g Spunlace Material	30*40cm	4			
		5	Standard Surgical Gown 标准型手术衣	40g SMS, Ultrasonic welding, knitted cuff, verco neck, book folding	XL/130*150cm	2			
		6	Suture bag 幘袋	80g Kraft paper	16*28cm	1			
		7	Op-tape背胶料	-	10*50	1			
		8	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190	1			
		9	Wrapper 外包布	30g SMS	100*100cm	1			
T. U. R. Pack 泌尿科手术包		1	TUR Drape 主单	40g SMS	170/240/175cm	1	1set/Blister bag, 10sets/Ctn, EO Sterile,Ctn size: 60x40x45cm	8.404USD/PC	2000
		2	Hand Towels 擦手纸	60g Spunlace Material,White	30*40cm	4			
		3	Standard Surgical Gown 标准型手术衣	40g SMS, Ultrasonic welding, knitted cuff, verco neck, book folding	XL/130*150cm	1			
		4	Reinforced Surgical Gown 加强手术衣	40g SMS, Ultrasonic welding, reinforced surgical gown, knitted cuff, verco neck, book folding, reinforced pad: 40g PP+PE film	XL/130*150cm	1			
		5	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190cm	1			
		6	Wrapper 外包布	30g SMS	100*100cm	1			
膝关节包/Knee Arthroscopy Pack-FJB0007		1	Knee Arthroscopy T-Drape 主单	40g SMS	200/250*315cm	1	1set/Head bag, 5bags/Ctn, EO Sterile, Ctn size: 60x31x43cm	14.54USD/PC	2000
		2	Split Drape/U型单	40g SMS, U hole, Adhesive	150*180cm	1			
		3	Stockinette 软质脚套	Elastic membrane + knitted liner	30*90cm	1			
		4	Mayo Stand Cover美亚袋	Width 80cm, 5.5s PE film+ reinforced 30g Hydrophilic PP 60*100cm	80*145cm	1			
		5	Hand Towels 擦手纸	60g Spunlace Material,White	30*40cm	2			
		6	Reinforced Surgical Gown 加强手术衣	40g SMS, Ultrasonic welding, reinforced surgical gown, knitted cuff, verco neck, book folding, reinforced pad: 40g PP+PE film	L/120*150cm	1			
		7	Hand Towels 擦手纸	60g Spunlace Material,White	30*40cm	2			
		8	Reinforced Surgical Gown 加强手术衣	40g SMS, Ultrasonic welding, reinforced surgical gown, knitted cuff, verco neck, book folding, reinforced pad: 40g PP+PE film	XL/130*150cm	1			
		9	Op-tape背胶料	-	10*50	1			
		10	Suture bag 幘袋	80g Kraft paper	16*28cm	1			
		11	Yankauer Suction Tube 杨克管	-	23+304cm	1			
		12	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190	1			
		13	Wrapper 外包布	30g SMS	100*100cm	1			
膝关节包/Knee Arthroscopy Pack-FJB0008		1	Knee Arthroscopy T-Drape 主单	40g SMS	200*300cm	1	1set/Head bag, 10sets/Ctn, Ctn size: 59*43*42cm	10.128USD/PC	2000
		2	Split Drape/U型单	40g SMS, U hole, Adhesive	150*180cm	1			
		3	Mayo Stand Cover美亚袋	Width 80cm, 5.5s PE film+ reinforced 30g Hydrophilic PP 60*100cm	80*145cm	1			
		4	Op-tape背胶料	-	10*50	1			
		5	Suture bag 幘袋	80g Kraft paper	16*28cm	1			
		6	Legging 腿套	40g SMS	45*60cm	1			
		7	Drape Sheet 铺单	Blue,40g SMS	140*190cm	1			
		8	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190cm	1			
		9	Wrapper 外包布	30g SMS	100*100cm	1			
Gynaecology ABD Set/妇科 ABD套装		1	Hand Towels 擦手纸	60g Spunlace Material,White	30*40cm	4	1set/blister bag 8bags/ctn Eo sterile 44*32*34cm	9.377USD/PC	2000
		2	Mayo Stand Cover美亚袋	Width 80cm, 5.5s PE film+ reinforced 30g Hydrophilic PP 60*100cm	80*145cm	1			
		3	Main Drape 主单	T Sharp+Medical Film+ Reinforced Pad+Adhesive	230*250cm	1			
		4	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190cm	1			

By-Pass
Pack 冠脉
搭桥手术
包



1	Cardiovascular U Split Drape/主单	40g SMS	218/260x340cm	1	1件/顶头袋8件/箱, 58*41*42	\$13.39	1000
2	Bilateral Split Drape/双侧分体单	40g SMS	200*200cm	1			
3	Top Drape/头单	40g SMS, 240cm side Adhesive 5*80cm, Reinforced 80g 30*90cm	150*240	1			
4	Mayo Stand Cover美亚袋	Width 80cm, 5.5s PE film+ reinforced 30g Hydrophilic PP 60*20cm	80*145cm	1			
5	Hand Towels 擦手纸	60g Spunlace material/60g 水刺布	30*40cm	2			
6	Shoe Cover鞋套	PP	40*20cm	2			
7	Shoe Cover	PE	40*20cm	2			
8	Op-tape背胶料	-	10*50	2			
9	Suture bag 幘袋	80g Kraft paper	16*28cm	2			
10	Side Drape 加强型边单	40g SMS, 90cm side Adhesive 5*88cm, Reinforced 80g 30*90cm	75*90	1			
11	Instrument Table Cover 桌布	Width 150cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*230cm	150*230cm	1			
12	Instrument Table Cover 桌布	Width 140cm, 55g PE case film+Reinforced 30g Hydrophilic PP 60*190cm	140*190cm	1			
13	Wrapper 外包布	30g SMS	100*100cm	1			

- 1.The standard product, Tolerance of material: ±1-2 gram, Tolerance of size: ±1.5-2cm
- 2.Price Validity: 7 days
- 3.USD foreign exchange currency is 1: 6.60 based on the day's current exchange rate
- 4.Payment Terms: 30% Deposit before Order, 70% balance before shipment.
- 5.Delivery: Within 25-35 days after receiving the deposit