


Product specification

1. Product Surgical gown, reinforced
2. Product Code
 1. Medium 11.34.01SN
 2. Large 11.34.02SN
 3. XLarge 11.34.03SN
 4. XXXLarge (X Large/Long) 11.34.04SN
 5. XXXLarge (X Large/Long) 11.34.05SN
3. Intended use. Suitable for procedures with medium to high fluid
4. Manufactured within clean room conditions
5. Material list
 1. Gown material 40gms AS SMS
 2. Reinforcement PE / tissue
 3. Cuffs cotton
6. Properties Liquid and impervious reinforced sleeve and chest panels and low linting
7. Dimensions
 1. Final fold dimensions 17 x 22 cm \pm 1 cm
8. Other features
 1. Velcro neck fastening
 2. Belt size 65cm
 3. Ultrasonic weld seams
 4. Tie-card showing shaded sterile / white non-sterile areas
 5. Coloured collar – White
9. Sterile packaging
 1. Handtowel 2 per gown 40x30 cm
 2. Wrap 20g SMS 60 x 60 cm
 3. Outer case 66*37*40cm 50gowns
 4. Sterilisation ETO
 5. Shelf life 3 years
10. Finished Product Dimensions



		11.34.01SN	11.34.02SN	11.34.03SN	11.34.04SN	11.34.05SN
		Medium	Large	XLarge	XXLarge (X Large/ Long)	XXXLarge (X Large/ Long)
A	Neck to hem ± 25 mm	1080	1120	1160	1200	1300
B	Chest ± 25 mm	650	660	670	700	710
C	Hem width ± 20 mm	1380	1420	1460	1500	1600
D	Sleeve length ± 20 mm	590	590	590	590	600
E	Cuff length ± 20 mm	70	70	70	70	70
F	Shoulder length ± 20 mm	200	200	210	210	220
G	Sleeve width ± 20 mm	280	290	290	290	300
H	Shoulder to hem ± 20 mm	1000	1040	1080	1120	1220
I	Neck width ± 20 mm	258	258	258	258	258
J	Sleeve reinforcement ± 20 mm	300	300	300	300	300
K	Chest reinforcement width ± 20 mm	500	500	500	500	500
L	Chest reinforcement width ± 20 mm	700	700	700	700	700