


Technical Data Sheet

Product specification

1. Product Patient Gown
2. Product Code Size F
FGPG-SMS35
3. Intended use. Suitable for procedures with low to medium fluid for Patient.
4. Manufactured within clean room conditions
5. Material list
 1. Gown material 35gms AS SMS
 2. Cuffs 3/4 Cuff
6. Properties Breathable, Antistatic, Liquid repellent, low linting
7. Dimensions
 1. Final fold dimensions 25*30 cm \pm 2 cm
8. Other features
 1. neck ties fastening
 2. Belt size 3.3*180cm
 3. Thread welding seams
 4. Color: Light Blue
9. Non Sterile packaging
 1. Inner PE Bag 10pcs/PE Bag
 2. Sterilization Non Sterile
 3. Shelf life 3 years
 4. Outer box 51*31*34cm/100pcs
10. Finished Product Dimensions



FGPG-SMS35

Free size(cm)

A	Shoulder width \pm 2cm	21
B	Chest \pm 2cm	68
C	Sleeve Length \pm 2cm	40
D	Axillary width	28
E	Neckline length	54
F	Distance from neck to Hem	105
G	Distance from shoulder to Hem	115
H	Bottom Width	137